

Easy use, optimal imaging.

Imaging Module IM 900®



Digital Slit
Lamp Imaging

 HAAG-STREIT
INTERNATIONAL

Tradition and Innovation



The intelligent choice for optimal imaging.


- Easy intuitive imaging
- Still imaging and high resolution video¹⁾
- Freeze Technology – high speed capture and data transfer
- Easy Touch – precise image control
- DFC – depth of field control
- History Trigger – never miss the perfect image
- Compact

¹⁾ Video capture is only available with EyeCap® SL

Imaging Module IM 900® the fully integrated compact imaging solution for the BQ 900® slit lamp.


Combination of the IM 900® with EyeCap™ or EyeSuite™ Imaging completes the perfect imaging system.

FREEZE TECHNOLOGY



Freeze technology: Captures at the precise moment you press the trigger.

HISTORY TRIGGER



Rewind time and freeze the moment: The history trigger not only captures one image but records the last few seconds in real time and lets you return to select the optimum image.

EASY TOUCH



Control the image: Easy touch more than just a trigger, allowing simple, ergonomic management of the camera settings.

DEPTH OF FIELD CONTROL (DFC)



True aperture: A unique feature of the Haag-Streit imaging systems, facilitating complete control of the depth of field.



From the leader in slit lamp imaging

Haag-Streit AG
Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone +41 31 978 01 11
Fax +41 31 978 02 82
info@haag-streit.com
www.haag-streit.com

 **HAAG-STREIT**
INTERNATIONAL

Tradition and Innovation