



**Prepare,  
Operate, Recover  
on one Chair.**  
MOBILE TABLE ROTATION

# **500 XLE** *surgiLine*

THE EFFICIENT OR EQUIPMENT

**SURGICAL CHAIRS & TABLES** FOR EYE- & MICROSURGERY

# 500 XLE

*The perfect reclining treatment stool for rotation systems in outpatient operating rooms.*



- Optimized price-performance ratio
- Ample legroom and unobstructed access to the patient's head
- Space-saving design for flexible use, outstanding maneuverability even in the most confined of spaces
- Multi-functional manual switch with 4 memory positions and auto-run function
- Automatic switch-off in the stand-by mode with reactivation button
- Sterile working conditions without interruption thanks to adjustment functions provided by the foot switches
- Not location-bound thanks to battery operation
- With equipotential bonding connection



- Easy maneuverability thanks to mechanic guidance operated via foot switch
- 360-degree rotating, electrically conductive double casters for outstanding movement characteristics
- Stool & table combination – no need for shifting or relocating the patients
- Service-friendly access for easy cleaning and maintenance
- Acoustic signal to warn against exhaustive battery discharge



- Easy mounting and dismantling, low access height
- Comfortable armrests for high patient safety, removable
- Comfortable upholstery in anthracite 197 standard shade



OPHTHALMIC  
SURGERY



multifunctional hand control



detachable foot control

## Standard equipment



*Power supply-independent battery box with integrated emergency shut-down switch.*



*Removable foot switch for infinite electronic height and head adjustment switch.*



*Standard sculptured headrest for high-precision patient positioning, adjustable via manual switches or foot switches.*



*2 mechanically operated brake levers with steering guide function.*



*Memory manual switch with four freely storable memory positions and auto-run function.*



*Telescopic stainless-steel accessories storage system.*

**SURGICAL CHAIRS & TABLES FOR EYE- & MICROSURGERY**

## **500 XLE** comfort

*The reclining treatment stool designed for a particularly high level of comfort. Our tried and tested designed in a whole new look.*



OPHTHALMIC  
SURGERY



- Two-tone comfort upholstery (anthracite 197/chrome standard shade) for even more comfortable sitting/reclining positions
- Standard color anthracite 197/chrome
- integrated knee support cushion



# surgiLine

*The modern surgeon's chair.*



## Armrests

- Rugged, 360-degree movable, and height-adjustable

## Seat System

- Seat system evenly distributes your weight as you make the slightest adjustment in position
- Anatomically formed, ergonomic backrest reduces pressure on inter vertebral discs and provides ideal support for the lumbar spine.

## Chassis

- 4 double casters
- Central brake acting on 4 wheels
- Power supply-independent battery box with integrated emergency shut-down switch
- Foot switch panel for moving the seat up & down



## Armrests

- Rugged, 360-degree movable, and height-adjustable

## Seat System

- Our unique dynamic seating system supports correct sitting postures while protecting your back

## Chassis

- Power supply-independent battery box with integrated emergency shut-down switch
- 4 double casters
- Central brake acting on 4 wheels
- 5th wheel for added resistance to tipping
- Foot switch panel for moving the seat up & down



# surgiTrend

*Designed with preventive medicine in mind.*



## *Designed for your health.*

As surgeons, the work you and your team have to perform is extremely demanding. Create optimal working conditions and avoid putting the entire surgical team under unnecessary stress by setting up a workplace that is individually adapted to suit your requirements.

# Your work is our challenge!

MOBILE TABLE ROTATION



**UFSK**  
INTERNATIONAL

**OSYS**

UFSK-INTERNATIONAL OSYS GmbH  
Kirchhoffstr. 1  
93055 Regensburg  
Germany

phone	+49 (0) 941 78862 – 15
fax	+49 (0) 941 78862 – 35
e-mail	<a href="mailto:info@ufsk-osys.com">info@ufsk-osys.com</a>
web	<a href="http://ufsk-osys.com">ufsk-osys.com</a>

**MADE IN GERMANY**

Design | Construction | Manufacturing