




Movies



Education



Crowding Bars




Contrast




*ProVideo*

Setting the STANDARDS for Vision Testing.



Pointer



ETDRS

**INNOVA**  
SYSTEMS, INC.

## PATTI PICS VALIDATED CHILDREN'S SYMBOLS

The Patti Pics optotype test has been designed to meet acuity test standardization criteria. They adhere to the Snellen principle and are calibrated to provide the same acuities as Sloan letters. Patti Pics are made up of five well-known symbols modified to provide a universal pediatric chart.



## AUTOLINK ADVANCED TESTING & INTEGRATION

INNOVA ProVideo® goes well beyond visual acuity to improve your practice – with advanced testing, exam room integration and superior optotypes. All of which translates into improved quality of care and efficiency in your practice.



## DISTANCE STEREO ADVANCED VISION TESTING

ProVideo STEREO® allows you to test stereopsis, establish accurate threshold and accurately monitor the patient's progress, making it useful for sports vision as well as pediatric and adult exams. ProVideo STEREO® provides excellent isolation using red green glasses. Stereo images include random dot, forced choice tests, and an animated butterfly.



## ProVideo CLASSIC INTEGRATED VISION TESTING

SOFTWARE FUNCTIONS + available as an upgrade item	ACUITY CHARTS
DataWall * Software Compatibility	Snellen, HOTV, Numbers & Allen Symbols
Customizable Chart Presentation Features	Tumbling E's / Landolt C Rings
VideoGate * Patient Education Interface	Validated ETDRS and Logmar Charts
DVD Interface	Worth 4 Dot
ProVideo * STEREO module +	Red / Green Screen
EMR Data Transfer +	Single Optotype / Single, Triple, Vertical Line
AutoLINK™ auto phoropter interface +	Crowding Bars @ 50%
Picture in Picture function	Fixation Dot
ADA Hearing Impaired Bars	Astigmatic Dials
NEI ATS & E-ETDRS Automated Test Protocols +	Optotype Pointer
Bi-Lingual Fixation Videos +	Size Display
	Cycling Lines/Optotypes
	Sloan Letters, Alternate Character Set
	White on Black Optotypes
	Fixation Disparity
<b>DISTANCE STEREO TESTS</b>	Patti Pics Validated Children's Symbols +
Automated Scripted Testing	Crowding Bars @ 100%
Animated Butterfly Pre-Verbal Test	Scaleable Custom Slides
Forced Choice Random Dots, Rings & Symbols	Contrast Flash Function
	Low Contrast Letters

**INNOVA**